

Screening protocol for the identification of spontaneous mutations in *Ciona intestinalis*

1. Acclimate individuals of wild *Ciona intestinalis* in 35-50 l tanks for 2-3 days (feeding regime as in Cirino et al. 2002): check daily animal conditions and discard dead/unhealthy individuals; if gonoducts contain few gametes, expose animals to continuous light for 2-3 days.
2. Dissect gonoducts and collect gametes separately in 2 ml vials:
 - open the animal at the base of the atrial siphon with a sterile blade;
 - collect eggs from the oviduct with a Pasteur pipette and transfer them in a 15 cm Petri dish filled with 0.22 μ M MFSW (millipore filtered natural sea-water);
 - collect sperm from the spermiduct with a Pasteur pipette and transfer it into a 2 ml vial in ice.
3. Perform *in vitro* self-fertilization by washing the pipette used for sperm collection into the egg-containing Petri dish (at 2-3 PM in order to have tadpole larvae ready for morphological screening the morning after) and store in an incubator at 17-18°.
4. Cryopreserve *Ciona* semen - see **Cryopreservation of *Ciona intestinalis* sperm**.
5. The morning after (9:30-10 AM), analyze carefully live progeny at the stereomicroscope for obvious phenotypes that affect gross anatomy. Stereotyped abnormal phenotypes must be around 20-25%. In case, add few menthol crystals for light anesthesia of tadpole larvae
6. If no obvious phenotypes are detected, fix about 100 larvae for ISH analysis - see **Automated whole-mount *in situ* hybridization on developmental stages of *Ciona intestinalis* for the identification of recessive mutations with subtle phenotype** -.
 - fix in 2 ml vials 1-1.5 hr in 1:1 8% PFA in 0.1 M MOPS pH 7.4:MFSW;
 - wash 3 x 10 min in 1 ml 1X PBT in DEPC H₂O;
 - dehydrate in EtOH series (3 x 10' in 1 ml 30%, 50% and 70% EtOH in DEPC H₂O);
 - change last 70% EtOH in DEPC H₂O and store at -20°C.

Fabio Crocetta, Paolo Sordino

Stazione Zoologica Anton Dohrn (SZN),
Naples, Italy

Reagents: 0.22 μ M Millipore filtered natural sea water (MFSW), DEPC - H₂O, paraformaldehyde, NaCl, MOPS, 10x PBS in DEPC-H₂O, EtOH, TWEEN-20

Equipment: blades (flame sterilized with alcohol), Pasteur pipettes, 15 cm Petri dishes, 2 ml vials, 50 ml Falcon tubes, 65 C water bath, Parafilm, 0.22 μ M Filters, Bunsen burner, incubator, 35-50 l tanks, stereomicroscope.

Additional information:

Cirino, P., A. Toscano, D. Caramiello, A. Macina, V. Miraglia, and A. Monte. Nov. 27, 2002. Laboratory culture of the ascidian *Ciona intestinalis* (L.): a model system for molecular developmental biology research. *Mar. Mod. Elec. Rec.* [serial online].
<http://www.mbl.edu/html/BB/MMER/CI/CirTit.html>

SOLUTIONS

8% PFA in 0.1 M MOPS pH 7.4

Reagent	Quantity (for 50 ml)	Final concentration
Paraformaldehyde	4 grams of Paraformaldehyde	
NaCl	5 ml of 5 M NaCl stock	0.5 M
MOPS	5 ml of 1 M MOPS pH 7.4 stock	0.1 M
DEPC H ₂ O	up to 50 ml	

Seal the 50 ml Falcon with Parafilm and mix well by shaking it manually.

Incubate in 65° water bath for 1 hr (shake manually each 10 minutes until paraformaldehyde is dissolved).

Filter the solution at 0.22 µM.

Aliquot in 2 ml vials and store at -20°.

1X PBT in DEPC H₂O

Reagent	Quantity (for 50 ml)	Final concentration
PBS	5 ml of 10X PBS stock	1X
100% TWEEN-20	50 µl	0.1%
DEPC H ₂ O	up to 50 ml	

Mix gently and store at room temperature.

30% ETOH in DEPC H₂O

Reagent	Quantity (for 50 ml)	Final concentration
100% EtOH	15 ml	30%
DEPC H ₂ O	up to 50 ml	

Mix gently and store at room temperature.

50% ETOH in DEPC H₂O

Reagent	Quantity (for 50 ml)	Final concentration
100% EtOH	25 ml	50%
DEPC H ₂ O	up to 50 ml	

Mix gently and store at room temperature.

70% ETOH in DEPC H₂O

Reagent	Quantity (for 50 ml)	Final concentration
100% EtOH	35 ml	70%
DEPC H ₂ O	up to 50 ml	

Mix gently and store at room temperature.
